

Reagents:

- His Dynabeads (or equivalent for your tag of interest)
- Yeast total RNA (10 mg/mL I have a lot of this)
- BSA (ultrapure, 50 mg/mL)
- 5M NaCl
- 1M Tris pH 8.
- 1M sodium phosphate pH 8
- Molecular biology grade (DNAse, RNAse, protease free) water.
- 4M imidazole pH 8 (or other concentrated eluent solution - glutathione, SDS, etc)

Buffers, made fresh from 5x stocks of PBST and TBST

PBST: 50 mM sodium phosphate pH 8.0, 300 mM NaCl, 0.01% Tween-20 (\pm 5 mM DTT) (make 1 mL/sample/selection)

TBST: 50 mM Tris-HCl pH 8, 150 mM NaCl, 0.05% Tween-20 (\pm 5 mM DTT) (1.1 mL/sample/selection)

Blocking: TBST with 0.1mg/mL BSA, 0.6 mg/mL yeast RNA 15 min @ 4°C. (150 μ L/sample/selection)

Library: 1 pmol/50 μ L, or 20 nM, in blocking buffer. (50 μ L/sample/selection)

Elution: PBST with 300 mM imidazole (50 μ L/sample/selection)

Selection:

1. Wash 25 μ L resin (volume for His dynabeads) with 2x 300 μ L PBST.
2. Dilute 5-25 μ g protein to 300 μ L in PBST. Incubate with resin 30 min @ 4°C. (Cold room rocker, tilted just short of vertical so the tubes are almost horizontal, set @ ~20%. Some may clump, just resuspend manually).
3. Save the protein FT that didn't bind, freeze on -80. Wash with 2x 200 μ L TBST.
4. Blocking: 100 μ L TBST with 0.1mg/mL BSA, 0.6 mg/mL yeast tRNA. 15 min @ 4°C.
5. Incubate 1-2 pmol library in 50 μ L TBST w/ BSA/tRNA, 60 min @ 4°C. Collect FT.
6. Wash 3x 200 μ L TBST. Switch tubes each for each wash (resuspend, take up suspension to new tube, place on magnetic rack and remove supernatant). This prevents contamination from library in/outside the tube.
7. Elute with 50 μ L PBST with 300 mM imidazole (5 min)

notes:

- store FTs/beads (samples with high library concentration) SEPARATE from the eluates. Freeze on dry ice.
- once protein is loaded onto beads, keep on ice as much as possible (until elution).